

MSH Emergency Department

MSH ED Tonopen Instructions

Last reviewed: 11-21-22

Locations:

1. Zone C Pyxis
2. Zone B – in eye cabinet in the Ophthalmology box
 - a. Code: **1710**

Instructions:



Instructional
Video



User Manual



iCare IC100 *the next gen tonopen*

1. Click any button to turn the tonometer on.
2. Take a new disposable probe, remove the cap and insert into the inferior probe base.
3. Have the patient sit up and stare ahead.
4. Place the device in the horizontal position with the superior pad touching the patient's forehead. The inferior probe should be 4-8mm away from the cornea.
5. Check that the lights surrounding the triangle button are green. If the color is red, then make sure your positioning is correct.
6. Click the triangle button, it will now make 6 separate measurements.
7. You will now hear a long beep and the measurement will be displayed on the screen.