

DOWNTIME REPORT VIEWING VIA PACS IN ZFP

1) Launch PACS in Chrome:

In-Hospital Network:

<https://mshzfp.rad.mountsinai.org/ZFP>

Out-of-Hospital Network:

<https://remotepacs.mountsinai.org/ZFP/>

GE Healthcare

Centricity™
Universal Viewer
Zero Footprint Client

User ID:
Password:

2) Enter **PACS** User/Pass (email # RADIOLOGY IT GROUP if needed as its different from MSHS email & Epic User/Pass)

Generic User ID: edzfp

Generic Password: edzfp1

Study Search REMOTEPACS

Filter

Patient Name: Last,First MI

Patient ID:

Accession #:

Date Of Birth: dd-MMM-yyyy

Referring Physician: Last,First MI

Exam Status:

Study Description:

Location:

Referring Services:

Modality: ☒ All

☐ CR ☐ DX ☐ CT ☐ PET/CT

☐ MR ☐ US ☐ XA

☐ PT ☐ NM ☐ GM

☐ MG ☐ OCT

☐ SC

Date: Yesterday Last 7 days

Search

Name	Patient ID	Description	Accession #
		SL CHEST 1 VIEW PORTABLE	
		SL CHEST 1 VIEW PORTABLE	
		SL FOOT LEFT	
		SL CHEST PA AND LEFT LA...	
		SL CHEST 1 VIEW PORTABLE	
		Unspecified CT ABDOMEN	
		SL CHEST 1 VIEW PORTABLE	
		SL CTA CHEST ABDOMEN P...	

3) Click report icon to launch report/case
Note it may be "unspecified"

Studies

Modality: All

2019

19 Nov CT CPACS Unspecified CT ABDOMEN

19 Nov CR CPACS SL CHEST 1 VIEW PORTA...

19 Nov CR CPACS SL CHEST 1 VIEW POR...

19 Nov CR CPACS SL CHEST 1 VIEW POR...

Series

Reports / Notes

CT Se:1 Topogram 0.6 T20f

CT Se:2 ART AX GATE 40%

CT Se:3

Unspecified CT ABDOMEN 19-Nov-2019

PRIMARY

test

4) Downtime Report in PACS

Se:1 Im: 1/1 H R L F

Se:2 Im: 1/85 A R L P

Se:3 Im: 1/85 A R L P

Se:4 Im: 1/85 A R L P